

Combat CARL offers all functions for TCCC and TECC training.


Life-size (183cm), 75 kg, with highly mobile extremities, no sensitive electronics, no gimmicks: Ideal for treatment and transport in extreme conditions. Suitable for practicing the standard algorithms.

## **MARCH**

- M Three included deep wounds can be applied to 7 sites and connected to bleeding modules. Wound packing is possible for all wounds, as is the application of tourniquets.
  - Airway management with oro- and nasopharyngeal as well as supraglottic devices.
- A Thoracic relief in monaldi and Bülau position

  Homostasis and access to the vascular system: i.v. in the force
- Hemostasis and access to the vascular system: i.v. in the forearm, i.o. sternal,
- R humerus head and proximal tibia with multiple-use pads. Voice output from prerecorded files or directly.
- **c** Robust and insensitive to water, dirt, cold and heat.

Н





Airway management nasopharyngeal



Airway management supraglottisch



Patient voice



Exchangeable eyes



Wounds in the thoracic region



Relief puncture according to Monaldi




Sternal intraosseous approach



Wound-Packing




i.V. additions



**Chest Drainage** 



IO access Humerus



IO access Tibia



TQ - placement



**Bandaging techniques** 



